



## FOOT & ANKLE PODIATRY INSTRUMENTS



Podiatry Instruments



Gelpi & Weitlaner  
Retractors



Hohmann Retractors



Orthopedic  
Instruments



Hip Surgery  
Instruments



Knee Surgery  
Instruments



Spine Instruments



Small Bone  
Instruments



Surgical Scissors



Surgical Elevators  
Phenolic Handle



Bone Holding  
Clamps



Needle Holders

Get in Touch  
**516 997 2180**  
[www.GerMedUSA.com](http://www.GerMedUSA.com)  
[sales@germedusa.com](mailto:sales@germedusa.com)

If we don't have  
the instrument you're looking for,  
we will **MANUFACTURE** it for you!

2417 Jericho Turnpike, #333  
Garden City Park, NY 11040 (USA)





## PODIATRY INSTRUMENTS

Podiatric or podiatry surgical instruments are suitable for removing callus, corns, nail trimming & cutting. These medical tools are ideal for structural deformities of the foot and ankle.

- Knife Instruments
- Nail Instruments
- Packer and Excavator
- Podiatry Bone Rongeurs

GerMedUSA is a trusted manufacturer and supplier of podiatry surgical instruments with multiple unique variations. Our high-quality instruments give precise results and enable the world's top care profession to manage the surgical industry by delivering quality care to patients.

## GELPI & WEITLANER RETRACTORS

The Gelpi Retractors are used for holding back tissue and organs during spine surgical procedures. The locking mechanism keeps the tool in place and provides exceptional hold for the surgery.

The Weitlaner Retractor with Solid Blades is designed to retract soft tissues, providing a clear view of the surgical site for the surgeon.

**These retractors have the following patterns:**

- Standard Gelpi Retractors
- Crossover Gelpi Retractor
- Deep Gelpi Lateral Retractor
- Gelpi Retractor With Speedlock
- Gelpi Retractor (Offset Arm Right)
- Weitlaner Retractor
- Solid Blades Retractor
- Gelpi Crossover Retractor (Neuroma) Deep Angle



## HOHMANN RETRACTORS

The Hohmann Retractors are also known as bony lever retractors. They are widely used in minimally invasive hip and knee surgeries.

**Have a look at the following patterns:**

- Bent Hohmann Retractor
- Mini Hohmann Retractors
- Hohmann Retractor Angled
- Right Angle Hohmann Retractor

The tip profile is leaf-shaped, round, or square for positive engagement with the bone. The contoured blade aids in easy exposure and insertion while minimizing the soft tissue damage.

## ORTHOPEDIC INSTRUMENTS

Orthopedic surgery instruments are designed for surgeries related to muscles, joints, ligaments, and bones. Our extensive range of orthopedic instruments accommodates a variety of medical approaches.

**Instruments used in orthopedic surgery are:**

- Bone Cutting Forceps
- Bone Files
- Bone Chisel
- Drill Guide
- Surgical Rongeurs
- Bone Rasps
- Femoral Ligament Cutter
- Periosteal Elevators
- Self Retaining Retractors
- Wire And Pin Management

Orthopedic surgical equipment is classified into different categories, and surgeons can pick the right instruments according to their choice.





## HIP SURGERY INSTRUMENTS

Hip Surgery Instruments are helpful for accurate bone placement and many other bone procedures such as fracture and bone manipulations.

These surgical instruments are crafted from premium-grade German surgical stainless material. The high-quality metal made our instruments durable and robust.

### Here are the types:

- Femoral Calcar Retractor
- Anterolateral Acetabular Rim Retractor
- Proximal Femoral Hook
- Greater Trochanteric Retractor

## KNEE SURGERY INSTRUMENTS

The knee Surgery Instruments are used to fix several knee joint problems. The surgical procedure may involve partial or complete knee replacement.

Our surgery instruments are the best high-precision products manufactured with great care in bulk quantities.

### Here are the types:

- Chandler Retractors
- Cobra Retractors
- Hohmann Retractors
- Knee Retractors
- Surgical Spreaders




## SPINE INSTRUMENTS

Spine surgical instruments are designed to treat spine surgeries. Our comprehensive spine and neurosurgery instrument portfolio includes hand-held and many self-retaining instruments made from German surgical stainless material. These spine instruments work efficiently and are considered the best surgeon's extended hand for better visibility and visualization.

## SMALL BONE INSTRUMENTS

Small bone procedures require precision, stability, and high-quality surgical tools to ensure successful outcomes. GerMedUSA offers a range of expertly crafted small bone instruments designed for optimal performance in orthopedic surgeries. Our instruments provide superior handling, durability, and efficiency, making them essential for delicate bone procedures.

- Bone Reduction Forceps
- Modified Mini Hohmann Retractors
- Pointed Fracture Reduction Clamp with Speed Lock
- Bone Fragment Clamps
- Weitlaner Retractor Sharp





## SURGICAL SCISSORS

Our multi-purpose surgical scissors are used by surgeons during an operation to cut tissues at the surface or inside the human body. Enrich your medical settings or clinics with our scissors that provide effective cutting-edge solutions. These surgical scissors are helpful for cutting.

Here are the types:

- Bandage Scissors
- Mayo Dissecting Scissors
- Metzenbaum Dissecting Scissors
- Microsurgery Scissors
- Operating Scissors
- Cardiovascular And Thoracic Scissors

## NEEDLE HOLDERS

Needle Holders are used during surgical operations to hold suture needles. The clamping mechanism provides a firm grip over the needle and locks it in place. They are similar in structure to a hemostat.

Here are the types:

- Mayo Hegar Needle Holder
- Mathieu Needle Holder
- Castroviejo Needle Holder
- DeBakey Needle Holder
- Derf Needle Holder
- Halsey Needle Holder
- Heaney Needle Holder
- Tungsten Carbide Needle Holder



## SURGICAL ELEVATORS PHENOLIC HANDLE

Periosteal elevators are specially designed to assist orthopedic surgeons in lifting and detaching the periosteum, the tough membrane that surrounds and protects bones. They are available in multiple designs and sizes, with different blade patterns for specific needs.

- Periosteal Elevators Round Edge
- Shoulder Penetrating Awl 7" Small Curved
- Joseph Rasp 6 3/4" Slight Curved 3mm
- 7 3/4" Wide Round Edge
- Cottle Elevator Joseph Curved

## BONE HOLDING CLAMPS

The German-forged bone clamps are specially designed to hold and stabilize bones during orthopedic procedures. We offer different sizes and designs to accommodate different types of bones and surgical techniques.

- Ikuta Bone Clamp
- Ludloff/Mau Osteotomy Fixation Clamp
- Fracture Reduction Clamp
- Bone Fragment Clamps Large
- Pointed Fracture Reduction Clamp

